



# SUPRAMALLEOLAR AFO

Patient \_\_\_\_\_ PO# \_\_\_\_\_ Diagnosis \_\_\_\_\_  
 Clinic \_\_\_\_\_ Due Date (Rush Fee May Apply) \_\_\_\_\_  
 Clinician \_\_\_\_\_ Height \_\_\_\_\_ Weight \_\_\_\_\_  
 Phone \_\_\_\_\_ Age \_\_\_\_\_ Sex \_\_\_\_\_ Shoe Size \_\_\_\_\_

<b>PATIENT ID</b>		<b>INTERNAL USE - DO NOT FILL IN PATIENT ID</b>
-------------------	--	---

## Step 1: Device Type

L B R <b>SMO with Plantarflexion Stop</b>	L B R <b>Standard SMO</b>
Alternate Footplate Plastic Lengths* <input type="checkbox"/> Proximal to Metatarsal Heads <input type="checkbox"/> Full Length   <input type="checkbox"/> EVA Footplate Padding (Full Length) *Standard Plastic Length is to Sulcus	Alternate Footplate Plastic Lengths* <input type="checkbox"/> Proximal to Metatarsal heads <input type="checkbox"/> Full Length   <input type="checkbox"/> EVA Footplate Padding (Full Length) *Standard Plastic Length is to Sulcus

Straps and pads included, but not shown above. Check back of form for material specifications. Plastic will be selected based on patient weight.

## Step 2: Modifications and Additions

Primary Color or Transfer Pattern:  Natural  Black  Blue  Pink  US Flag  Forest Camo  Other

ANKLE AND REARFOOT	ANKLE AND REARFOOT	MIDFOOT	FOREFOOT
L B R <b>Do Not Change Ankle*</b> Or Change to _____	L B R <b>Medial Wedge</b> Height or Angle _____	L B R <b>Increase Arch Height</b>	L B R <b>Do Not Change*</b> Or Change to _____
L B R <b>Do Not Change Rearfoot*</b> Or Change to _____	L B R <b>Lateral Wedge</b> Height or Angle _____	L B R <b>Decrease Arch Height</b>	L B R <b>Metatarsal Pad</b> Small Medium Large
	L B R <b>Posterior Wedge</b> Height or Angle _____	L B R <b>Navicular Relief</b>	L B R <b>Toe Crest Pad</b>
		L B R <b>Cuboid Relief</b>	

\*Ankle, rearfoot, and forefoot correction to neutral is standard.

## Step 3: Order Notes and Special Requests




RECEIVED DATE	COMPLETED DATE

# SUPRAMALLEOLAR AFO

## Step 4: Measurements

Measurements are not necessary, but will improve the accuracy of the final product.



## Material Specifications

Call with questions or for other material options.  
Warranty information available at [www.AODCFAB.com](http://www.AODCFAB.com)

### SMO with Plantarflexion Stop

- Polypropylene Homopolymer Shell and Posting
- 35 Durometer EVA Padded Footplate *optional*
- 25 Durometer EVA Malleolar and Arch Padding
- 35 Durometer EVA Posterior Stop Padding

### SMO

- Polypropylene Homopolymer Shell and Posting
- 35 Durometer EVA Padded Footplate *optional*
- 25 Durometer EVA Malleolar and Arch Padding

INTERNAL USE - DO NOT WRITE BELOW

CUSTOMER COMMUNICATIONS